



University of
Salford
MANCHESTER

Ram locus is a key regulator to trigger multidrug resistance in *Enterobacter aerogenes*

Molitor, A, James, CE, Fanning, S, Pages, JM and Davin-Regli, A

<http://dx.doi.org/10.1099/jmm.0.000667>

Title	Ram locus is a key regulator to trigger multidrug resistance in <i>Enterobacter aerogenes</i>
Authors	Molitor, A, James, CE, Fanning, S, Pages, JM and Davin-Regli, A
Publication title	Journal of Medical Microbiology
Publisher	Microbiology Society
Type	Article
USIR URL	This version is available at: http://usir.salford.ac.uk/id/eprint/44909/
Published Date	2018

USIR is a digital collection of the research output of the University of Salford. Where copyright permits, full text material held in the repository is made freely available online and can be read, downloaded and copied for non-commercial private study or research purposes. Please check the manuscript for any further copyright restrictions.

For more information, including our policy and submission procedure, please contact the Repository Team at: library-research@salford.ac.uk.

Target region	Primer	Sequence	Ann Temp
For sequencing			
<i>marR</i>	M5	CCT-GTT-TAT-TAC-GCT-CGG-CGT	64°C
	M4IIB	CCA-GCG-ACA-ATG-GAG-ATT-CCA-GG	64°C
<i>marA</i>	M7	TAT-GAT-TGA-AAT-CAA-ACG-GCG	64°C
	M12	CAG-ACG-AAG-TGG-TCA-TGC-TTG	64°C
<i>ramR</i>	RamR1	TGA-AAG-GGT-AGG-TTA-GCG	64°C
	RamR2	CAG-ACG-AAA-GCC-CCA-TCC	64°C
<i>ramA</i>	Rafud	GGT-ACC-GTT-TAG-TCA-CGC-GAA-TG	64°C
	RamARevI	CCG-CTT-GGA-AAC-GCT-GTT-ATC	64°C

Supplementary Table 1: Primers for sequencing

Name	Sequence		Annealing Temp	
<i>omp35</i>	5'	GGCGGCGACGGCTGGA ACTACAC	3'	65°C
	5'	GCAAAGCTCAGGCCGTCCACCAGA		
<i>omp36</i>	5'	AGACCTACAACGCAACTCGCACTG	3'	60°C
	5'	AACGGCCCATGTCTTTACCTTTAG		
<i>acrA</i>	5'	ATCTACCGCGCTGGCTACCGTTCA	3'	61°C
	5'	TCGCCTTACCGTTTTCTTGCTTCA		
<i>ramA</i>	5'	ACATTCGCGAGCGTAAACTGC	3'	58°C
	5'	CCCGGCGGCTGATTGAAGGTA		
<i>ramR</i>	5'	AGATGCGCTTAGTACAGACAATG	3'	53°C
	5'	AACCGCTTTCACCGACAG		
<i>marA</i>	5'	AAGTTTCGGAGCGCTCAGGTTAC	3'	58°C
	5'	AGTTTTTGC GCAATCTCCGTCAG		
<i>marR</i>	5'	CCAATTTCCGGGTGCTCTGTTCC	3'	57°C
	5'	TGCCGCGCTTGTCATTGG		
<i>soxR</i>	5'	GCTGCTGTGCGCTGGGGAAGTG	3'	62°C
	5'	GTCGCCAGCGGAATACCAATACG		

Supplementary Table 2 : Primers for quantitative realtime PCR

Target region	Primer	Sequence	Ann Temp
For cloning			
<i>ramR</i>	RamR(Fw)XhoI	GTT-ACG-CTC-GAG-GTG-GCT-CGT-CCT-AAG-AGT	64°C
	RamR(Rev)SacI	GTT-ACG-GAG-CTC-TTA-CGA-CTC-ATT-CTC-ATG	64°C
	RamR(Rev)SacI EA27	GTT-ACG-GAG-CTC-TTA-CGA-CTC-ACT-CTC-ATG	64°C
<i>ramA</i>	RamA(Fw)BamHI	CCT-GTA-GGA-TCC-CAT-TTA-ACG-CCT-GGT-GGC-GC	64°C
	RamA(Rev)EcoRI-2	GTT-ACG-GAA-TTC-TCA-GTG-CGC-GCG-GCT-GTG	64°C
For control			
	T7	TAA-TAC-GAC-TCA-CTA-TAG-GG	64°C
	SP6	ATT-TAG-GTG-ACA-CTA-TAG	64°C

Supplementary Table 3: Primers for cloning experiments