



University of  
**Salford**  
MANCHESTER

# Ram locus is a key regulator to trigger multidrug resistance in *Enterobacter aerogenes*

Molitor, A, James, CE, Fanning, S, Pages, JM and Davin-Regli, A

<http://dx.doi.org/10.1099/jmm.0.000667>

<b>Title</b>	Ram locus is a key regulator to trigger multidrug resistance in <i>Enterobacter aerogenes</i>
<b>Authors</b>	Molitor, A, James, CE, Fanning, S, Pages, JM and Davin-Regli, A
<b>Type</b>	Article
<b>URL</b>	This version is available at: <a href="http://usir.salford.ac.uk/id/eprint/44909/">http://usir.salford.ac.uk/id/eprint/44909/</a>
<b>Published Date</b>	2018

USIR is a digital collection of the research output of the University of Salford. Where copyright permits, full text material held in the repository is made freely available online and can be read, downloaded and copied for non-commercial private study or research purposes. Please check the manuscript for any further copyright restrictions.

For more information, including our policy and submission procedure, please contact the Repository Team at: [usir@salford.ac.uk](mailto:usir@salford.ac.uk).

<b>Target region</b>	<b>Primer</b>	<b>Sequence</b>	<b>Ann Temp</b>
<b>For sequencing</b>			
<i>marR</i>	M5	CCT-GTT-TAT-TAC-GCT-CGG-CGT	64°C
	M4IIB	CCA-GCG-ACA-ATG-GAG-ATT-CCA-GG	64°C
<i>marA</i>	M7	TAT-GAT-TGA-AAT-CAA-ACG-GCG	64°C
	M12	CAG-ACG-AAG-TGG-TCA-TGC-TTG	64°C
<i>ramR</i>	RamR1	TGA-AAG-GGT-AGG-TTA-GCG	64°C
	RamR2	CAG-ACG-AAA-GCC-CCA-TCC	64°C
<i>ramA</i>	Rafud	GGT-ACC-GTT-TAG-TCA-CGC-GAA-TG	64°C
	RamARevI	CCG-CTT-GGA-AAC-GCT-GTT-ATC	64°C

**Supplementary Table 1:** Primers for sequencing

<b>Name</b>	<b>Sequence</b>		<b>Annealing Temp</b>	
<i>omp35</i>	5'	GGCGGCGACGGCTGGAACACTACAC	3'	65°C
	5'	GCAAAGCTCAGGCCGTCCACCAGA		
<i>omp36</i>	5'	AGACCTACAACGCAACTCGCACTG	3'	60°C
	5'	AACGGCCCATGTCTTTACCTTTAG		
<i>acrA</i>	5'	ATCTACCGCGCTGGCTACCGTTCA	3'	61°C
	5'	TCGCCTTACCGTTTTCTTGCTTCA		
<i>ramA</i>	5'	ACATTCGCGAGCGTAAACTGC	3'	58°C
	5'	CCCGGCGGCTGATTGAAGGTA		
<i>ramR</i>	5'	AGATGCGCTTAGTACAGACAATG	3'	53°C
	5'	AACCGCTTTCACCGACAG		
<i>marA</i>	5'	AAGTTTCGGAGCGCTCAGGTTAC	3'	58°C
	5'	AGTTTTTGCGCAATCTCCGTCAG		
<i>marR</i>	5'	CCAATTTCTGGGTGCTCTGTTCC	3'	57°C
	5'	TGCCGCGCTTGTCATTGG		
<i>soxR</i>	5'	GCTGCTGTGCGCTGGGGAAGTG	3'	62°C
	5'	GTCGCCAGCGGAATACCAATACG		

**Supplementary Table 2 :** Primers for quantitative realtime PCR

<b>Target region</b>	<b>Primer</b>	<b>Sequence</b>	<b>Ann Temp</b>
<b>For cloning</b>			
<i>ramR</i>	RamR(Fw)XhoI	GTT-ACG-CTC-GAG-GTG-GCT-CGT-CCT-AAG-AGT	64°C
	RamR(Rev)SacI	GTT-ACG-GAG-CTC-TTA-CGA-CTC-ATT-CTC-ATG	64°C
	RamR(Rev)SacI EA27	GTT-ACG-GAG-CTC-TTA-CGA-CTC-ACT-CTC-ATG	64°C
<i>ramA</i>	RamA(Fw)BamHI	CCT-GTA-GGA-TCC-CAT-TTA-ACG-CCT-GGT-GGC-GC	64°C
	RamA(Rev)EcoRI-2	GTT-ACG-GAA-TTC-TCA-GTG-CGC-GCG-GCT-GTG	64°C
<b>For control</b>			
	T7	TAA-TAC-GAC-TCA-CTA-TAG-GG	64°C
	SP6	ATT-TAG-GTG-ACA-CTA-TAG	64°C

**Supplementary Table 3:** Primers for cloning experiments